

Monoclonal Anti-Eribulin Antibody, Rabbit IgG (1M1G11)

Catalog # ERN-MY2012b



BIOSYSTEMS
Acro

Source

Monoclonal Anti-Eribulin Antibody, Rabbit IgG (1M1G11) is a Rabbit monoclonal antibody recombinantly expressed from HEK293 cells.

Clone

1M1G11

Species

Rabbit

Isotype

Rabbit IgG | Rabbit Kappa

Conjugate

Unconjugated

Antibody Type

Recombinant Monoclonal

Reactivity

Chemical

Immunogen

Eribulin.

Specificity

Specifically recognizes Eurlbin.

Application

Application	Recommended Usage
ELISA	0.06-31 ng/mL

Purity

>90% as determined by SDS-PAGE.

Purification

Protein A purified/ Protein G purified

Formulation

Lyophilized from 0.22 μ m filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

SDS-PAGE



Discounts, Gifts,
and more!



➔ www.acrobiosystems.com

9/10/2024

Monoclonal Anti-Eribulin Antibody, Rabbit IgG (1M1G11)

Catalog # ERN-MY2012b




BIOSYSTEMS
Acro


Monoclonal Anti-Eribulin Antibody, Rabbit IgG (1M1G11) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).

Bioactivity-ELISA

Monoclonal Anti-Eribulin Antibody, Rabbit IgG (1M1G11)-Bridging ELISA
0.1 µg of Monoclonal Anti-Eribulin Antibody, Rabbit IgG (1M1G11) per well



Monoclonal Anti-Eribulin Antibody, Rabbit IgG (1M1G11) ELISA
0.1 µg of Farletuzumab Ecteribulin per well



Immobilized Monoclonal Anti-Eribulin Antibody, Rabbit IgG (1M1G11) (Cat. No. ERN-MY2012b) at 1 µg/mL (100 µL/well) can bind Farletuzumab Ecteribulin diluted with 100% Mouse Serum with a linear range of 0.02-1.25 µg/mL (QC tested).

Immobilized Farletuzumab Ecteribulin at 1 µg/mL (100 µL/well) can bind Monoclonal Anti-Eribulin Antibody, Rabbit IgG (1M1G11) (Cat. No. ERN-MY2012b) with a linear range of 0.06-2 ng/mL (Routinely tested).

Background

Eribulin is a synthetic analogue of the macrocyclic polyether halichondrin B, which was originally isolated from the Asian sea sponge *Halichondria okadai*. Eribulin binds specifically to the β -tubulin subunit on the (+) end of the microtubule and potently inhibits elongation of the formed microtubule, while having little or no effect on microtubule depolymerization. Eribulin's potent antimitotic activity and nonmitotic effects on tumor biology make it an interesting candidate for investigation as a MTA payload for ADCs.

Clinical and Translational Updates

Discounts, Gifts,
and more!



www.acrobiosystems.com

9/10/2024