

PE-Labeled Human CD7 Protein, His Tag (Site-specific conjugation)

Catalog # CD7-HP2E3



Synonym

CD7,GP40,TP41,LEU-9,TP40

Source

PE-Labeled Human CD7, His Tag (CD7-HP2E3) is produced via site-specific conjugation of PE to Human CD7 (26-180), His Tag under optimal conditions with a proprietary technology. Human CD7, His Tag (CD7-HP2E3) is expressed from human 293 cells (HEK293). It contains AA Ala 26 - Pro 180 (Accession # [P09564-1](#)).

Predicted N-terminus: Ala 26

Molecular Characterization

CD7(Ala 26 - Pro 180)
P09564-1 Poly-his

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 20.0 kDa.

Conjugate

PE

Excitation Wavelength: 488 nm / 561 nm

Emission Wavelength: 575 nm

Application

Please note that this product is NOT compatible to streptavidin detection system.

Formulation

Lyophilized from 0.22 μ m filtered solution in PBS, 0.5% BSA, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage


For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please protect from light and avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.

Bioactivity-FACS



5×10^5 of anti-CD7 CAR-293 cells were stained with 100 μ L of 1:50 dilution (2 μ L stock solution in 100 μ L FACS buffer) of PE-Labeled Human CD7 Protein, His Tag (Cat. No. CD7-HP2E3) and negative control protein respectively. PE signal was used to evaluate the binding activity (QC tested).

Discounts, Gifts,
and more!



➤ www.acrobiosystems.com

12/29/2023

PE-Labeled Human CD7 Protein, His Tag (Site-specific conjugation)

Catalog # CD7-HP2E3



BIOSYSTEMS
Acro

Background

T-cell antigen CD7 (CD7) is also known as GP40, LEU-9, TP41 and Tp40. CD7 is a protein that in humans is encoded by the CD7 gene, this gene encodes a transmembrane protein which is a member of the immunoglobulin superfamily. CD7 has been shown to interact with PIK3R1. This protein is found on thymocytes and mature T cells. It plays an essential role in T-cell interactions and also in T-cell/B-cell interaction during early lymphoid development.

Clinical and Translational Updates

Please contact us via TechSupport@acrobiosystems.com if you have any question on this product.

Discounts, Gifts,
and more!



➤ www.acrobiosystems.com

12/29/2023