

Monoclonal Anti-HSV-2-Glycoprotein D Antibody, Human IgG1 (3H9)

Catalog # GLD-M609



Source

Monoclonal Anti-HSV-2-Glycoprotein D Antibody, Human IgG1 (3H9) is a chimeric monoclonal antibody recombinantly expressed from HEK293, which combines the variable region of a mouse monoclonal antibody with Human constant domain.

Clone

3H9

Isotype

Human IgG1 | Human Kappa

Conjugate

Unconjugated

Antibody Type

Recombinant Monoclonal

Reactivity

Human

Immunogen

Recombinant HSV-2 (strain 333) Envelope Glycoprotein D (gD) derived from human 293 cells.

Specificity

This product is a specific antibody specifically reacts with Glycoprotein D of HSV-2 and HSV-1.

Application

Application	Recommended Usage
ELISA	1-500 ng/mL

Purity

>90% as determined by SDS-PAGE.

Purification

Protein A purified/ Protein G purified

Formulation

Lyophilized from 0.22 μ m filtered solution in PBS, pH7.4 with trehalose as protectant.

Contact us for customized product form or formulation.

Reconstitution

Please see Certificate of Analysis for specific instructions.

For best performance, we strongly recommend you to follow the reconstitution protocol provided in the CoA.

Storage

For long term storage, the product should be stored at lyophilized state at -20°C or lower.

Please avoid repeated freeze-thaw cycles.

This product is stable after storage at:

- -20°C to -70°C for 12 months in lyophilized state;
- -70°C for 3 months under sterile conditions after reconstitution.


SDS-PAGE

Discounts, Gifts,
and more!



Monoclonal Anti-HSV-2-Glycoprotein D Antibody, Human IgG1 (3H9)


Catalog # GLD-M609



Monoclonal Anti-HSV-2-Glycoprotein D Antibody, Human IgG1 (3H9) on SDS-PAGE under reducing (R) condition. The gel was stained with Coomassie Blue. The purity of the protein is greater than 90% (With [Star Ribbon Pre-stained Protein Marker](#)).


Bioactivity-ELISA

Monoclonal Anti-HSV-2-Glycoprotein D Antibody, Human IgG1 (3H9) ELISA
0.1 µg of HSV-1 (strain 17) Envelope Glycoprotein D (gD), His Tag per well



Monoclonal Anti-HSV-2-Glycoprotein D Antibody, Human IgG1 (3H9) Conc. (ng/mL)

Monoclonal Anti-HSV-2-Glycoprotein D Antibody, Human IgG1 (3H9) ELISA
0.1 µg of HSV-2 (strain HG52) Envelope Glycoprotein D (gD), His Tag per well



Monoclonal Anti-HSV-2-Glycoprotein D Antibody, Human IgG1 (3H9) Conc. (ng/mL)

Immobilized HSV-1 (strain 17) Envelope Glycoprotein D (gD), His Tag (Cat. No. GLD-V52H3) at 1 µg/mL (100 µL/well) can bind Monoclonal Anti-HSV-2-Glycoprotein D Antibody, Human IgG1 (3H9) (Cat. No. GLD-M609) with a linear range of 1-31 ng/mL (QC tested).

Immobilized HSV-2 (strain HG52) Envelope Glycoprotein D (gD), His Tag (Cat. No. GLD-V52H4) at 1 µg/mL (100 µL/well) can bind Monoclonal Anti-HSV-2-Glycoprotein D Antibody, Human IgG1 (3H9) (Cat. No. GLD-M609) with a linear range of 0.1-2 ng/mL (Routinely tested).

Background

Herpesvirus infections are widely spread throughout the world population. Herpes simplex virus (HSV) belongs to the α -herpesvirus subfamily. There are two main types of HSV, HSV-1 and HSV-2, which infect humans. HSV-2 mainly causes genital lesions, whereas HSV-1 is involved in both oral and genital infections. Glycoprotein D (gD) is a structural component of the herpes simplex virus type 1 (HSV-1) envelope which is essential for virus entry and fusion with host cells. gD plays an important role by binding to the host receptors such as herpes virus entry mediator (HVEM) and nectin-1, a member of the immunoglobulin (Ig)-like cell adhesion molecules.

Clinical and Translational Updates

Discounts, Gifts,
and more!

